Chord-X® Mitral Chordal Replacement System Instructions for Use

EN	English Instructions for Use.....	2
BG	Български Инструкции за употреба.....	9
CS	Čeština Návod k použití.....	16
DA	Dansk Brugsanvisning.....	23
DE	Deutsch Gebrauchsanweisung.....	30
EL	Ελληνικά Οδηγίες χρήσης.....	37
ES	Español Instrucciones de uso.....	44
ET	Eestikeelne kasutusjuhend.....	51
FA	فارسی دستور العمل استفاده.....	58
FI	Suomi Käyttöohjeet.....	65
FR	Français Mode d'emploi.....	72
HR	Hrvatski – Upute za uporabu.....	79
HU	Magyar Használati utasítás.....	86
ID	Bahasa Indonesia Instruksi Penggunaan.....	93
IT	Italiano Istruzioni per l'uso.....	100
JA	日本語 取扱説明書.....	107
KO	한국어 사용설명서.....	114
LT	Lietuvių k. Naudojimo instrukcijos.....	121
LV	Lietošanas norādījumi latviešu valodā.....	128
NL	Nederlands Gebruiksaanwijzing.....	135
NO	Norsk Bruksanvisning.....	142
PL	Polski Instrukcja użycia.....	149
PT	Português Instruções de utilização.....	156
PT-BR	Português do Brasil Instruções de uso.....	163
RO	Română Instrucţiuni de utilizare în limba.....	170
RU	Русский Инструкция по применению.....	177
SL	Slovenska navodila za uporabo.....	184
SR	Srpski Uputstvo za upotrebu.....	191
SV	Svensk bruksanvisning.....	198
TR	Türkçe Kullanım Talimatları.....	205
UK	Українська Інструкція для застосування.....	212
VI	Tiếng Việt Hướng Dẫn Sử Dụng.....	219
ZH	中文 使用说明书.....	226

Artivion, Inc.
1655 Roberts Blvd. NW
Kennesaw, GA 30144, USA
+1-770-419-3355
www.artivion.com

Jotec GmbH
Lotzenäcker
72379 Hechingen
Germany
+49 (0) 7471 922 0
www.jotec.com

OMCS Medical GmbH
Aegeristrasse 5
6300 Zug
Switzerland

OMC Medical Limited
Planet House,
North Heath Lane
Horsham, West Sussex
RH12 5QE
United Kingdom

Australian Sponsor
Device Technologies
1 Garigal Rd.
Belrose NSW 2085 Australia
+61 2 9972 8203

1. OVERVIEW OF THE CHORD-X MITRAL CHORDAL REPLACEMENT SYSTEM

The Chord-X® Mitral Chordal Replacement System (Chord-X System) is a premade artificial suture system comprised of the Chord-X® Suture System or Chord-X® Pre-Measured Loops, and Chord-X® Chordal Sizer. The Chord-X System is intended for use by trained cardiac surgeons with experience in performing mitral valve repair procedures.

1.1 Intended Purpose

The Chord-X Suture System or Chord-X Pre-Measured Loops are intended for the replacement of chordae tendinae during a mitral valve repair procedure. The device is intended to be used in patients with mitral regurgitation due to defective native mitral valve chords (i.e., ruptured or elongated native chords). The Chordal Sizer is intended to assist the surgeon in determining the length of the chordal tendinae to be replaced, for selection of the Chord-X Loop size (pre-measured or adjusted).

1.2 DEVICE DESCRIPTION

The Chord-X Pre-Measured Loops (Figure 1) and the Chord-X Suture System (Figure 2) are comprised of nonabsorbable, monofilament 2-0 suture strands manufactured from expanded polytetrafluoroethylene (ePTFE) which have a total length of approximately 32 inches (81cm). The sutures are attached to a pledget made of polytetrafluoroethylene (PTFE). Each suture strand is threaded through an AISI 302 stainless steel curved needle in either 3/8 circle (135° arc) or 1/2 circle (180° arc) point-to-point distance, also referred to as 'circle' size, with an 18mm taper point (for a total of 8 needles). The Chord-X System includes three (3) chordal suture pairs either with pre-measured loops and pre-tied knots (Chord-X Pre-Measured Loops) or untied, adjustable suture strands (Chord-X Suture System) with attached needles for suturing to the mitral leaflets, one (1) pair of suture strands with an attached pledget and needles for the papillary muscle attachment, and one (1) extra pledget. Each suture pair is positioned on a series of baseplates for suture management and protection during transport. If fewer chordal sutures are desired, excess segments can be removed without jeopardizing the rest of the prosthesis. The Chord-X Pre-Measured Loops and Chord-X Suture System are in a box containing five (5) single-packaged cases (with a Patient Implant Card included in each case) which are provided sterile by ethylene oxide gas for single use only.


Figure 1. Chord-X Pre-Measured Loops


Figure 2. Chord-X Suture System

The Chord-X Chordal Sizer (Figure 3) is designed to assist surgeons in determining the proper length of the chordal tendineae to be replaced for selection of the Chord-X Pre-Measured Loop size (or adjustment of the loop size in the Chord-X Suture System). Use of the Chord-X Sizer is optional. There are multiple chordae tendineae measuring techniques and instruments and surgeons typically have their own preference (Nasso et al. 2023). The Chord-X Chordal Sizer is constructed of Polycarbonate/Acrylonitrile Butadiene Styrene (PC/ABS) thermoplastic with a polyester label. Five (5) Chord-X Chordal Sizers are provided within each box. Each sizer is for single use only and is supplied STERILE within a double-pouched system. Each Chord-X Chordal Sizer is sterilized by gamma irradiation.


Figure 3. Chordal Sizer

The available configurations of the Chord-X System are provided in Table 1 below.

Table 1: Chord-X Mitral Chordal Replacement System

Chord-X System	Loop Length	3/8 Circle Needle	1/2 Circle Needle
Chord-X Suture System	Adjustable	CXL-20-1838-0	CXL-20-1812-0
Chord-X Pre-Measured Loops	12mm	CXL-20-1838-12	CXL-20-1812-12
	16mm	CXL-20-1838-16	CXL-20-1812-16
	20mm	CXL-20-1838-20	CXL-20-1812-20
	24mm	CXL-20-1838-24	CXL-20-1812-24
Chord-X Chordal Sizer	N/A	CXCS	

The performance characteristics for the Chord-X Suture System and Chord-X Pre-measured loops are defined by United States <USP29> and European <EP9.0> pharmacopoeia requirements for Nonabsorbable Surgical Sterile Sutures, to which the Chord-X System is compliant. The performance characteristics of the Chord-X Chordal Sizer include locking force to extend the chordal sizer.

- The locking force or force to extend the pointer is ≥ 0.5 lbs and ≤ 4 lbs.
- Rebound of the pointer after removing thumb pressure is less than 0.5mm.
- The scale is accurate to ± 1 mm.

The qualitative and quantitative information on the implant materials is listed below:

Table 2: Chord-X Materials

Chord-X Component	Material	Mass	Mass Composition (%)
Suture	Expanded Polytetrafluoroethylene "ePTFE" thread	0.30g	79%
Pledgets	Polytetrafluoroethylene "PTFE"	0.08g	21%

The ePTFE Suture has been expanded to produce a porous microstructure which is approximately 50% air by volume. Additionally, the suture is undyed and contains no additives. The needles supplied for attaching the ePTFE chords contain cobalt (CAS No. 7440-48-4; EC No. 2311580), which is defined as CMR 1B in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

2. INDICATIONS FOR USE

The Chord-X System is indicated for treatment of patients with mitral regurgitation who are not at high surgical risk as judged by a multidisciplinary heart team.

3. CONTRAINDICATIONS

The Chord-X System is contraindicated for use in ophthalmic surgery, microsurgery, and peripheral neural tissue. Additionally, the Chord-X System is contraindicated in patients with acute severe MR due to papillary muscle rupture who are recommended valve replacement in accordance with current AHA/ACC and ESC/EACTS guidelines, and patients with infectious endocarditis (as endocarditis can weaken the leaflet tissue and affect the stability of the suture).

4. WARNINGS

Physicians should inform patients about all potential risks and adverse events discussed in these instructions for use.

- **Do not** implant the prosthesis if the device is visibly damaged or altered, or the expiration (use-by) date has expired, as this can adversely affect the sterility of the device.
- **Use:** The Chord-X Pre-Measured Loops and Chord-X Suture System should only be implanted by physicians who have received appropriate training and are experienced in cardiac surgical techniques, including mitral valve repair.
- **Accurate Sizing Required:** The replacement of chordae tendineae requires accurate sizing to ensure satisfactory functioning of the valve mechanism post repair or replacement. Chordae tendineae replaced with sutures that are either too long or too short can result in failure to correct valvular incompetence resulting in extended bypass time, surgical delay, regurgitation, heart failure, atrial fibrillation, stroke, and pulmonary hypertension. (See Directions for Use in **Section 9** for the recommended sizing techniques to determine correct length.)
- **Do Not Re-sterilize:** The prosthesis has been sterilized using ethylene oxide. Do not re-sterilize.
- **Do Not Reuse:** The Chord-X Pre-Measured Loops and the Chord-X Suture System prostheses are designed for single use only; do not reuse the device. Reuse may cause device failure or procedural complications including device damage, compromised device biocompatibility, and device contamination. Reuse may result in infection, serious injury, or patient death.
- **Potential Tissue In-Growth:** In rare case, the use of the prosthesis may result in suture calcification and/or the attachment of the suture to the tissue it penetrates which may result in difficult removal of the sutures and/or regurgitation/leakage of the mitral valve.

5. PRECAUTIONS

Caution should be taken in treatment of patients with mitral regurgitation due to rheumatic valve disease, extensive valve prolapse, and leaflet or annular calcification, for which the benefits of complex repair vs. mitral valve replacement should be weighed.

- Administration of anticoagulation therapy should be individualized for each patient per institutional guidelines, to avoid potential thromboembolic or bleeding events.
- Mishandling of this suture prosthesis, like any other suture, can result in severe injury or death to the patient. Care should be taken when handling the device in order to avoid any inadvertent damage to the suture material or pledgets, which may impact its longevity (See Directions for Use in **Section 9** for additional information). Additionally, care should be taken when implanting the device in order to avoid any inadvertent damage to native cardiac tissue.
- The Chord-X needles associated with the suture prosthesis are sharp. Special care should be taken when handling these needles to avoid inadvertent needlesticks. Additionally, used needles should be discarded with appropriate sharps containers and/or in accordance with the implant facility's policy for such medical wastes.
- When implanting the suture prosthesis, be aware that if two (2) strands are cut, there will be a free-floating piece of suture that will need to be retrieved.
- Shortening of the artificial chordae while knotting could result in overcorrection of the prolapsed leaflet and restrict the leaflet motion.
- Do not use metal knot pushers or knot pushers with imperfections such as burrs or flashing as it can damage the suture.

6. POTENTIAL ADVERSE EVENTS & UNDESIRABLE SIDE EFFECTS

Potential adverse events that may be associated with any mitral valve repair procedure include, but are not limited to:

- Allergic reaction
- Tissue reaction
- Tissue damage
- Prosthesis damage
- Recurring valve regurgitation/leakage
- New onset atrial fibrillation
- Thromboembolism
- Stroke
- Heart failure
- Permanent disability
- Reoperation for recurrent severe MR and/or explant of device
- Death

Potential adverse event that may be associated with use of the Chord-X System, but are not limited to:

- Thromboembolism (due to degradation of the pledget material)
- Rupture of the artificial chords
- Detachment of the chords from the papillary muscle or leaflet tissue, potentially leading to leaflet prolapse resulting in recurrent mitral regurgitation and/or decreased cardiac output

7. MAGNETIC RESONANCE IMAGING (MRI) SAFETY AND COMPATIBILITY

The ePTFE Sutures and Pledgets (the implantable portion of the Chord-X Pre-Measured Loops and Chord-X Suture System) are considered to be MR Safe.

8. PATIENT INFORMATION AND COUNSELING

Patients who have mitral valve dysfunction (e.g., moderate to severe mitral valve regurgitation) who require repair or replacement of their chordae tendineae, and who can tolerate a surgical heart valve repair procedure, should discuss their options for treatment using the Chord-X Mitral Chordal Replacement System with their physician. Patients undergoing mitral valve repair should be monitored in accordance with the ACC/AHA and/or ESC guidelines (e.g., biannually or as recurrent MR symptoms arise) after their initial procedure to ensure mitral valve function. The physicians should inform patients about all potential risks and adverse events discussed in the IFU.

9. DIRECTIONS FOR USE

- 9.1** Before opening the sterile package, examine it carefully for evidence of damage or compromised sterility. If intact, remove the device from the sterile packaging by locating the separated corner of the pouch (i.e., top right-hand corner) and pulling the two (2) layers in opposite directions. For Chordal Sizer Inner Pouch, locate the separated edge of the packaging and aseptically pull the two (2) layers in opposite directions to provide access to the device. Care should be taken to avoid damaging the device during use (e.g., bending the chordal sizer shaft or hooks). Aseptic technique should be practiced when handling the device during the implant procedure.

WARNING: DO NOT use if package or device is damaged. Damage to the sterile barrier may result in potential infection or fever.

- 9.2** As with any suture, care should be taken to avoid damage when handling.
- Avoid crushing or crimping the suture with surgical instruments or exposing the suture to sharp edges, which may adversely affect the tensile strength of the suture.
 - Avoid using knot pushers with imperfections such as burrs or flashing which may damage the suture.
 - Avoid untying the pre-tied suture knots on each pre-measured loop.
- 9.3** Unpack Chord-X Pre-Measured Loops or the Chord-X Suture System from the baseplate.
- After determining the appropriate loop length, choose the corresponding size of Chord-X Pre-Measured Loops or the Chord-X Suture System. Once selected, open the package using sterile aseptic technique and note that the prosthesis is positioned on a series of baseplates for suture management. Once opened, perform the following: Remove foam baseplates from the case.
 - Unwrap the papillary suture with attached pledget from the baseplate.

- Unwrap chordal prosthesis sutures 1, 2, and 3 from their respective baseplates.
- Remove the additional PTFE pledget from case.

Note: Both the Chord-X Pre-Measured Loops and the Chord-X Suture System prosthesis feature one (1) suture pair for a papillary muscle and three (3) independent pre-measured chordal suture loops or three (3) suture strands (for the leaflet) complete with needles. If fewer chordal sutures or strands are desired, excess segments are easily removed without jeopardizing the rest of the prosthesis.

CAUTION: Be aware if two (2) strands are cut, there will be a free-floating piece of suture that will need to be retrieved.

9.4 Determine the Chord-X Pre-Measured Loops or Chord-X Suture System loop length. Using the Chord-X Chordal Sizer, or other surgeon preferred technique, determine the length of the loop needed in order to maintain the coapting surfaces of the two leaflets at the same level in the ventricle. The following are recommended sizing technique with Chord-X Chordal Sizer.

- Use the blue leaflet hook to support the edge of the leaflet in the planned position.
- Use the thumb-slide to extend the papillary hook to the planned site of attachment to the papillary muscle (see Figure 4).
- The leaflet hook may be rotated to any desired angle relative to the papillary hook to improve visibility. Rotate the leaflet hook by using the grips at the base of the barrel.
- The colors of the scale match the label colors for the corresponding Chord-X Pre-Measured Loops prostheses.
- The scale number indicates the size of chordal loops produced by tying suture around the two hooks. The scale is accurate to ± 1 mm. Rotate the barrel so that the leaflet hook is aligned with the papillary hook prior to tying loops. After tying, pull the papillary hook back to remove tied loops from the hooks.
- Once the artificial chordae tendineae are in place, verify that the sizer remains intact, a saline test is also recommended to verify that the chordal loop size is correct and that the valve is competent.


Figure 4. Sizing

9.5 Remove extra loops (optional, as needed at surgeon discretion):

Remove undesired prosthetic loops by cutting a single strand of suture between the knot and the papillary suture strand and discard the undesired loop by removing the single strand in its entirety.

Figure 5. Removal of loops from Chord-X Pre-Measured Loops

Remove undesired suture strand by cutting a single chordal suture strand and discard the undesired loop by removing the single strand in its entirety.

Figure 6. Removal of Loops from Chord-X Suture System

9.6 Papillary Suture Placement: Place the needle of the papillary suture through the papillary muscle, with a pledget placed on each side of the muscle. One (1) pledget is pre-attached to the prosthesis and the second pledget is provided in the case. The papillary suture is then tied using at least seven (7) flat square throws; additional throws may be needed per physician discretion (Figure 7). The Needles are removed and discarded upon securing the ePTFE chords to the papillary muscle.



Figure 7. Papillary Suture Placement

9.7 Leaflet placement and knot tying:

Once the suture has been anchored to the papillary muscle, pass each of the two (2) strands of the Chord-X loop through the free edge of the prolapsed leaflet. Tie the sutures using at least seven (7) flat square throws, taking care not to entrap the native chordae as the suture is positioned and tied. Continue placement of the remaining Chord-X Pre-Measured Loops as needed to complete artificial chordal replacement, tying each one as described above.

Figure 8. Chord-X Pre-Measured Loops Leaflet Placement and Knot Tying

Once the suture has been anchored to the papillary muscle, pass the required chordal suture strands through the free edge of the prolapsing leaflet, taking care to not entrap normal chordae as the suture system is individually positioned and tied. Establish the length of each chordae tendineae replacement suture using the Chord-X Chordal Sizer, or other surgeon preferred method, so that the suture is set to maintain the coapting surfaces of the two leaflets at the same level in the ventricle.

Figure 9. Chord-X Suture System Leaflet Placement and Knot Tying

9.8 Surgeons may employ a variety of intraoperative tests to assess chordal length before closing the left atrium. In one such test, the ventricle is filled with fluid under pressure by insufflating it with saline, giving antegrade cardioplegia, or performing both maneuvers. With such testing, if residual regurgitant jets are identified, the surgeon may adjust the length of the strands.

9.9 Knot the Chord-X Suture System prosthesis suture ends using at least seven (7) flat square throws. Continue placement of the remaining Chord-X Suture System sutures strands as needed to complete artificial chordal replacement, tying each one as described above.

CAUTION: At least seven (7), equally tensioned, flat square throws are required to produce a secure knot when tying the Chord-X Pre-Measured Loops or the Chord-X Suture System for Mitral Chordal Replacement made of ePTFE suture. Additional throws may be necessary depending on surgical circumstances. When tying knots with the Chord-X Pre-Measured Loops or the Chord-X Suture System prostheses made of ePTFE suture, tension should be applied by pulling each strand of the suture in opposite directions with equal force. As the knot is tensioned, the air in the suture is forced out. Care should be taken to avoid using a jerking motion which could break the suture. Uneven tensioning of a well-formed square knot may result in an insecure knot, potentially impacting the effectiveness of the device. When the Chord-X Pre-Measured Loops or the Chord-X Suture System for Mitral Chordal Replacement is properly tensioned and formed, standard surgical knotting techniques will produce a secure knot.

CAUTION: Shortening of the artificial chordae while knotting could result in overcorrection of the prolapsed leaflet and restrict the leaflet motion.

9.10 Utilize surgeon preferred method of imaging to confirm proper loop length and proper coaptation of leaflets and successful valve repair. After the patient is off cardiopulmonary bypass, transesophageal echocardiography should be used to confirm the results of the hydrodynamic testing and evaluate the valve repair. Figure 10 depicts the finished Chordal Replacement.

Figure 10. Final Chordal Tendinae Replacement

10. CLINICAL BENEFIT

The clinical benefits of the Chord-X System include high survival rates with low incidence of cardiac, neurological, and visceral complications, and low recurrent severe Mitral Regurgitation with low rates of device failure and reoperation.

The summary of safety and clinical performance (SSCP) is available in the European database on medical devices (EUDAMED), where it is linked to the Chord-X Chordal Mitral Replacement System Basic UDI-DI (81357002CXLVB). The EUDAMED URL is <https://ec.europa.eu/tools/eudamed>.

11. REPORTING OF SERIOUS INCIDENTS AND ADVERSE EVENTS

Serious Incidents that occur in relation to the Chord-X procedure or device use should be reported to Artivion Field Assurance at 1-800-438-8285 or by email to fieldassurance@artivion.com, and to the national regulatory authority as applicable.

12. STORAGE, RETURNS, AND DISPOSAL

The sterility expiration date of the Chord-X devices is recorded on the outer package label. Appropriate inventory control should be maintained so that devices with earlier expiration dates are preferentially implanted and expiration is avoided. To protect the device, it should be stored in its outer box until used. Once used, the device should be discarded in accordance with the facilities policy for medical waste. The storage environment should be clean, cool, and dry.

Used needles should be discarded within appropriate sharps containers and/or in accordance with the implant facility's policy for such medical wastes. Explanted or unused portions of the prosthesis should be discarded in accordance with the implant facility's policy for medical waste.

13. PATIENT IMPLANT CARD

The Chord-X Pre-Measured Loops or Chord-X Suture System is packaged with an implant card, all patients should be instructed to keep this card in their possession at all times. Instructions for completion of Patient Implant Card (to be completed by the healthcare institution/provider):

1. Patient name (first, middle, last) or patient ID
2. Date of implantation (YYYY-MM-DD)
3. Name and address of the implanting healthcare institution/provider
4. Physician's name

The Patient Implant Card can also be printed from www.artivion.com/implantcard.

The patients should be instructed to refer to the patient information leaflet located at www.artivion.com/information-leaflet. The web address is located on the patient implant card. The explanation of all symbols is provided at the same weblink. Healthcare providers should print a copy of the patient information leaflet and explanation of symbols for patients who do not have easy access to Internet.

The current revisions of all On-X product IFUs may be found at: <https://www.artivion.com/eIFU>.

A printed version of the Instructions for Use will be provided within seven days upon request through any of the contact means listed below:

Phone: +1-770-419-3355 • Fax: 770-590-3753 • customerservice.us@artivion.com

14. SYMBOLS

	Batch Code
	Catalogue Number
	Caution
	CE Mark
	Consult Instructions for Use
	Contains Hazardous Substances – Cobalt
	Date of Implantation
	Date of Manufacturer and Country of Manufacturer
	Distributor
	Do Not Discard
	Do Not Re-sterilize
	Do Not Reuse
	Do Not Use if Package is Damaged
	Double Armed Needle
	Double Sterile Barrier System, Sterilized Using Gamma Irradiation (Chordal Sizer)
	European Authorized Representative
	Health Care Center or Doctor

	Importer
	Information Website for Patients
	Manufacturer
	Medical Device
	MR Safe
	Patient Identification
	Quantity
	Single Sterile Barrier with Protective Packaging, Sterilized Using Ethylene Oxide (Chord-X Loops and Sutures)
	Store in Cool, Dry Place
	Swiss Authorized Representative
	Taper Point Needle
	Translation
	U.S. federal law restricts this device to sale by or on the order of a physician.
	UK Authorize Representative
	UK Conformity Assessment
	Unique Device Identifier
	Use By Date

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